



NEAR VISION TEST CARD

	DISTANCE CORRELA- TION	JAEGER	PT	VISUAL EFF%
D T 4	20/800		72	5%
L E S 3	20/400		42	10%
R F X B N	20/250	18	30	15%
P O 5 7 A	20/200	16	26	20%
9 V M C L	20/100	10	14	50%
K S 3 Z 7	20/70	7	10	65%
N R E T X	20/50	5	8	75%
O R D F M P	20/40	3	6	85%
V J F X G H	20/30	2	5	90%
P 3 5 E A R	20/20	1	4	100%

This card has been specially designed for the vision care practitioner to aid in standardized measurements of near point acuity. This card should be held at approximately 16 inches away from the patients face under standard room illumination.